

Inclusion of attention gates in a 3D CNN model improves the performance in new T2 lesion quantification

M.Guirao¹, N.Nerseyan³, S.Valverde¹, R.Bramon¹, A.Clèrigues¹, L.Valencia¹, G.Álvarez³, M. Puig³, A. Gifreu³, A.Oliver², X.Lladó², À.Rovira⁴, L.Ramió-Torrentà³

¹ Tensor Medical, Girona, Spain.

² Department of Computer Architecture and Technology, University of Girona, Girona, Spain.

³ Neurodegeneration and Neuroinflammation Research Group, Girona Biomedical Research Institute (IDIBGI). Girona Multiple Sclerosis and Neuroimmunology Unit. Neurology Department. Dr. Josep Trueta University Hospital and Santa Caterina Hospital. Department of Medical Sciences. University of Girona

⁴ Section of Neuroradiology and MR Unit Department of Radiology, Vall d'Hebron Hospital Universitari, Barcelona, Spain.



tensormedical



Multiple Sclerosis



Every 5 minutes a new case of MS is diagnosed worldwide



MS affects more than **2.8M** worldwide



3 out of 4 of new diagnosed cases will be **women**



Annual global costs exceed **€ 145.000 M**

Multiple Sclerosis (Spain)



55.000 patients
diagnosed in Spain



Prevalence of **120**
patients / 100.000
habitants



A new case is
diagnosed every
5 hours



Annual global costs
exceed
€ 2300 M
(**€ 43.000 PPPY**)


Multiple Sclerosis (treatment)




Disease modifying
treatments (DMT)
€ 40.000 PPPY



DMT can slow down
the **progression** of
the disease



MRI biomarkers
evaluate the
response to
treatments



MRI brain lesions
are nowadays the
best **biomarker**

New / enlarging T2 lesions



Nowadays most neuradiologists perform **visual inspection**



Prone to **errors** due to **repositioning**, diffuse and **confluent** chronic lesions




Visual inspection decreases the **sensitivity** detecting lesion activity



Sub-optimal administration of the best **personalized** treatment

3D ResU-Net with Attention Gates (AGResU-Net)



3D Convolutional Neural Network

3D architecture

↑ Sensitivity

↑ Accuracy

Attention gates

Focus on relevant regions

Ignore irrelevant zones

MRI data

Acquired at **Josep Trueta** hospital (Girona).

100 subjects with 2 MRI (mean interval of 14 months [range 6 – 36 months]).


2 x **Philips Achieva 1.5T** scanner.

Four sequences using the **same protocol**:

- 3D-FLAIR
- 3D-T1-w
- T2-w
- PD-w

SAGITAL 3D-T1_MPRAGE
(TR/TE/TI=2300/2.98/900ms,
voxel size=1x1x1mm)


3D-FLAIR (TR/TE=5000/394ms,
voxel size=1x1x1mm)




AXIAL PD_T2_tse (TR/TE1/TE2=2500/16/91ms, voxel size=0.8x0.8x3mm)

Methodology: automated methods


Two different methods are used to determine the number of **new / enlarging** T2 lesions



M1 - Four input sequences



M2 - FLAIR only + Attention Gates



Methodology: evaluation

All methods are compared to a **gold standard** (reference outcome of new / enlarging lesion) that is composed as follows:



Gold standard

=



Visual
radiological
report

+



Automated
method


+



One expert radiologist
(> 15 years of experience)


Results: visual standard report vs. automated methods

Visual standard report (VSR) obtained less than 35% of sensitivity while M1 and M2 reached 76% and 90%, respectively



Results: visual standard report vs. automated methods

VSR detected **26 lesions** over 100 patients. **M1** raised the value to **56** and **M2** extended the count to **69**




Results: radiological activity

VSR identified **14 patients** with **26 NEDA-compatible lesions**

M1 identified **19 patients** with **59 lesions**

M2 found **24 patients** with **69 lesions**



Conclusions



Attention gates improved **detection** of T2 lesions in MS patients **significantly**



AGResU-Net with FLAIR sequences only **outperformed other methods**



Suitable to assist the radiologist **during the treatment** response in MS patients

Thanks

mguirao@tensormedical.ai



tensormedical

